

# Immunoassay test report



|          |         |          |        |                   |                  |
|----------|---------|----------|--------|-------------------|------------------|
| Patient: | ROVER   | Species: | Canine | Patient ID:       | 106182           |
| Client:  | CABIGAS | Gender:  | Male   | Sample No.:       | 01               |
| Doctor:  |         | Age:     |        | Time of analysis: | 2025/12/03 19:02 |

| Lab item   | Current result | Ref. Ranges  |
|------------|----------------|--------------|
| cNT-proBNP | 414            | pmol/L 0-900 |

Operator:

## Report Explan.

### cNT-proBNP

#### Result indications:

<900 pmol/L Normal

900-1800 pmol/L Suspected

>1800 pmol/L Abnormal

#### Clinical significance:

cNT-proBNP is cleaved from BNP which is produced by canine cardiomyocytes and increases with excessive stretching of the cells. cNT-proBNP concentration reflects the severity of cardiac activation after stimulation (e.g. stretching) and is the optimal biomarker for diagnosis and rule out heart failure. Used for prognostic evaluation of heart failure and therapeutic efficacy monitoring. Increase indicates an increased risk of heart failure.

Note: Due to the complexity and individuality of disease diagnosis, the report interpretation is only for your reference. Please consult your doctors for clinical diagnosis results.  
The results only applies to this test sample.

Test Instrument: Mindray vetXpert I3

Time of Printing: 2025-12-03 19:03:30



BATINGA ANIMAL MEDICAL CENTER  
SM CITY CDO UPTOWN BRANCH  
0906 121 1260 – 088 327 1837

Global Pioneer of Comprehensive Animal Medical Solutions  
Better healthcare for all - Since 1991

