

Immunoassay test report



Patient:	HATCHI	Species:	Canine	Patient ID:	119858
Client:	JALA	Gender:	Spayed	Sample No.:	01
Doctor:		Age:	10Y	Time of analysis:	2025/12/06 17:38

Lab item	Current result	Ref. Ranges	2025/11/06
cNT-proBNP	225	pmol/L 0-900	<div><div></div></div> 5914

Operator:

Report Explan.

cNT-proBNP

Result indications:

<900 pmol/L Normal

900-1800 pmol/L Suspected

>1800 pmol/L Abnormal

Clinical significance:

cNT-proBNP is cleaved from BNP which is produced by canine cardiomyocytes and increases with excessive stretching of the cells. cNT-proBNP concentration reflects the severity of cardiac activation after stimulation (e.g. stretching) and is the optimal biomarker for diagnosis and rule out heart failure. Used for prognostic evaluation of heart failure and therapeutic efficacy monitoring. Increase indicates an increased risk of heart failure.

Note: Due to the complexity and individuality of disease diagnosis, the report interpretation is only for your reference. Please consult your doctors for clinical diagnosis results.
The results only applies to this test sample.

Test Instrument: Mindray vetXpert I3 Time of Printing: 2025-12-06 17:39:11



BATINGA ANIMAL MEDICAL CENTER
SM CITY CDO UPTOWN BRANCH
0906 121 1260 – 088 327 1837

Global Pioneer of Comprehensive Animal Medical Solutions
Better healthcare for all - Since 1991

